

# BT-200 Fetal Doppler

# HI-bebe

**HI-bebe is showing you what fetus is saying.**

HI-bebe is pocket size doppler for Obstetricians, midwives and pregnant women.

- Compact and light
- High-Sensitivity Doppler probe
- Ergonomic design
- Large-Display & clean sound
- Easy to use
- Long time continuous use (6 hours)



Sound version



# BT-200 Fetal Doppler **Hi-bebe**

## Specifications

<b>Display</b>	• LCD Back light • Singnal quality display	• Low battery display
<b>Heart rate range</b>	• 50~240 bpm	
<b>FHR Accuracy</b>	• $\pm 2\%$	
<b>Ultrasound Frequency</b>	• 2 MHz	
<b>Ultrasound Intensity</b>	• $< 10 \text{ m W/cm}^2$	
<b>Speaker output power</b>	• 1.2 W (Built in Loudspeaker)	
<b>Audio Output for Earphone</b>		
<b>Auto Shut off Sound</b>	• 1 Minute	
<b>Auto Shut off Power</b>	• 5 Minutes	
<b>Power</b>	• 1.5 V Battery X 2 (Type : AA)	
<b>Sensitivity</b>	• 10-12 Weeks onward	
<b>PC Interface</b>	• Sound card (BCM-220)	
<b>Battery Life</b>	• 360 Minutes	
<b>Accessories</b>	• Carrying case • Gel (60ml)	• Battery (1.5V AA X 2)

## Waterproof (Optional)

## PC Software

### BCM220

- Record and play fetal heart sound (.wav format)

